

# ZC3

get back, faster



# FAST RECOVERY THERAPY UNIT

## COLD - HOT , AIR OR LIQUID COMPRESSION



### NO ICE / NO WATER

Therapeutic cooling technology without ice or water inside.

### CLOSED & CLEAN SYSTEM

ZAMAR DRIPLESS SYSTEM preloaded with glycol doesn't need to refill or drain off liquid.

### INNOVATIVE COOLING

Advance technology based on real-time temperature parameters.

### STAY FOCUSED ON TEMPERATURE

Set exact temperature you want to reach.

### INTERMITTENT COMPRESSION

Improves benefits of cold therapy and healing with 3 levels of air-compression or liquid.


### NEW SYSTEM PULSATION or COMPRESSION

Possibility to select AIR compression or liquid pulsation.

### INSPIRED BY YOUR SHAPE

Complete range of Anatomic wraps designed on your body.

## FLEXIBILITY & POWER



ZC3



# FLEXIBILITY & POWER



## WHEN EVOLUTION NEVER ENDS

### PLUG 'N' PLAY RECOVERY

No set up time and preparation needed.

### FRIENDLY USE

All treatment parameters manageable in one screen.

### REPEAT YOUR FAVORITE

Possibility to store up to 12 personal treatments for frequent use

### TWO STAGES OF TEMPERATURE

Ability to set 2 different temperatures and time phases



# ZT-CUBE V.3

# ZC3



## Features

<b>APPLICATION</b>	COLD & HOT
<b>TEMPERATURE RANGE</b>	+0°C +40°C (32°F 104°F)
<b>SCREEN</b>	LCD 7" color TouchScreen
<b>WRAPS</b>	Liquid / Air Wraps
<b>Glycol</b>	NON TOX freeze
<b>Air Compression levels</b>	mmHg: 50 - 75 - 100
<b>setup time limit</b>	300 h.
<b>Medical class for Wraps</b>	BF
<b>Number of Wraps</b>	up to 2
<b>Salvable therapies</b>	1 + 12
<b>DIMENSIONS</b>	cm. 37 x 32 x h33 inc. 15 x 12 x h13
<b>WEIGHT</b>	15 kg. - (33 lb.)

- LCD 7 inc, TFT color touch screen
- Easy software, Intuitive Menu icons
- Possibility to select two different temperature and treatment time
- Possibility to select two different types of massage, Air or Liquid
- Set programs that doctors can customize and save up to 12 memory areas
- Totally automated cycle management
- Fast connexion/detachment of wraps and pipes with Dripless connectors
- Anatomic Wraps pre charge with glycol
- NON-TOX cooling liquid



**ZAMAR MEDICAL d.o.o.**

Sv. Martin, 6 - 52450 - Vrsar - (HR)  
Tel: +385 (0)52 496 111 - VAT: 96642432252  
info@zamar.care

[www.zamar.care](http://www.zamar.care)

1370

ISO 13485  
BUREAU VERITAS  
Certification



ISO 9001  
BUREAU VERITAS  
Certification

